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GB Instructions for Use for Prosthetic Components 4plus6Line of implant system OT-F² for occlusal screw fixations

Introduction

Basic information and explanations for this implant system (such as implant diameter, color coding, etc.) are contained in the product catalog OT-F²/OT-F.

The surgical procedure for implant insertion and general indications and contra-indications are contained in the instructions for use for surgery of the implant systems.

The abutments of the 4plus6Line as a multi-unit-abutment serve for the restoration of edentulous jaws. Bar restorations are also possible with this abutment line. Please see also the indications for the ProFileLine abutments.

Illustrations, indications and article numbers to the prosthetic components mentioned below are also contained in detail in the product catalog OT-F²/OT-F3.

Malocclusion or overloading of the implant abutment connections, the abutments and the complete prosthetic construction should be strictly avoided. The pre-surgical planning of an implant restoration should include important aspects such as:

Implant positions, implant diameters, number of implants and the planned prosthetic design of the superstructure in addition to the bite situation of the patient. The distal extension should not exceed the width of a premolar.

The patient should be instructed to observe any changes in fixation of his denture during daily hygiene and to consult the dentist immediately, if required. We recommend a recall at regular intervals.

Important notes:

- The 4plus6Line abutments (0°/17°/39° in gingival heights 1.5 and 3.00 mm) are supplied in sterile condition, as they can be inserted into the inserted implants immediately. The color-coded final screw is mounted on the abutment. If the sterility is not guaranteed by incidental damage of the package, the preparation instructions for medical devices should be observed.
- These abutments contain an internal thread connection which avoids the loss of the abutment screw during surgery.
- The abutments should be fixed with 35 Ncm. The corresponding prosthetic driver are contained in the product catalog at the section "Drivers" (octagon driver for the one-piece non-antiroation 0° abutments and the 1.30 hex driver for the antiroation safe angulated abutments). The abutment screws must not be cemented or glued within the implant. The abutments should not be modified.
- All other prosthetic components are not supplied in sterile condition (gingival healing abutments, impression posts, final abutments, etc.). The package box of the healing abutments and the impression copings contains the Instructions for Preparation of Medical Devices.
- The prosthetic final abutments (cylinder POM, cylinder titanium) are supplied with the laboratory screw mounted (M1.4) and in addition with a brown coded final screw (M1.4).
- In principle, only new abutment screws which have not been used in the lab, may be used for final restoration in the mouth of the patient.
- Detailed information is contained in the product catalog in section "Screws".
- The implants of 3.80 mm diameter have an internal thread connection of M1.6, the implants of 4.10/5.00 mm of an internal thread of M1.8. Please observe this when ordering new components.
- The final screws are color-anodized for an improved identification. Please use only the final screw for final fixation of the abutments in the mouth of the patient.
- Color-coding according to the color guide system is contained in the 4plus6Line abutments in the region of the FourByFour connection (yellow 3.80, red 4.10, blue 5.00 mm diameter).
- Abutments for the implant diameter 3.40 mm are not available for this restoration concept.

The treatment concept / Indications

The treatment concept favors the rehabilitation of edentulous jaws and avoids

comprehensive augmentation procedures. Only four implants are required (in the maxilla possibly six), whereby the posterior implants of preferably 16 mm length should be inserted in approx. 30° angle towards distal. We recommend to use a drilling guide for the pilot bone preparations. This treatment method allows an immediate loading of the implants with a screw-on temporary bridge within only one visit.

The implants are provided immediately after insertion with the corresponding sterile 4plus6Line abutments, and after suturing an open impression is taken on abutment level. The correct selection of the angulation and gingival height and the direction of the occlusal screw-on fixation (rotation position of the implant) are of importance. The corresponding implant analogs are mounted and a plaster model is prepared. The temporary bridge prepared earlier should then be placed. At the beginning, only one titanium cylinder should be placed, and then the next following should be polymerized intraorally without tension. After finishing, the denture is placed on the fixed 4plus6Line abutments and screwed occlusally (brown color-coded M1.4 screw, 25 Ncm).

During the time between implant insertion and insertion of the temporary or final restoration, the 4plus6Line abutments should be protected by gingival healing abutments (PEEK).

Contraindications:

- Primary connections of abutments
- direct Cast-on technique
- Single tooth restoration
- Cemented bridge restorations

Important advice:

These abbreviated instructions do not replace a comprehensive explanation of the complete treatment concept. The user should be an experienced implantologist who has intensively dealt with the philosophy of this procedure and the special surgical and prosthetic protocol. A close and coordinated cooperation between dentist and dental technician is of decisive significance for the success of this treatment concept.

The following prosthetic components are available for the OT-F² implant system:

Sterile packed 4plus6Line abutments in the gingival heights of 1.50 and 3.00 mm for:

- Ø3.80 mm in angulations of 0° and 17°
- Ø4.10 mm in angulations 0°, 17° and 30°

Packed in unsterile condition:

- Gingival healing abutment [PEEK]
- Impression coping, open
- Implant analog
- Cylinder [Titanium]
- Cylinder [POM]

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OT medical GmbH