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Instructions for Prosthetic Components of Implant Systems OT-F² and OT-F³

Introduction

Illustrations, indications and product catalog numbers referring to the prosthetic components mentioned are contained in detail in the OT-F²/OT-F³ product catalog. Malocclusion and overloading of the implant abutment connection, the abutments and the complete prosthetic construction should be absolutely avoided.

Important instructions:

- All prosthetic components (except for the 4plus6Line) of the OT-F²/OT-F³ implant systems are supplied in non-sterile condition.
- Basically only new and non-used abutment screws (i.e. not used in the dental laboratory) should be used for final restoration in the mouth of the patient. The corresponding prosthetic abutments are supplied with a lab screw mounted and also a color coded final screw in addition.
- The implants of Ø 3.40/ Ø 3.80 mm diameter have an internal thread of size M1.6, the implants of Ø 4.10/ Ø 5.00 mm have an internal thread of M1.8. Please consider this for processing or re-ordering of the components.
- The final screws are color-coded for distinction. Please use only the final screw for final fixation of the abutments in the mouth of the patient.
- All abutments (except for CreativeLine and GoldLine) have an internal screw thread which prevents a possible loss of the abutment screw during processing.
- All abutments should be fixed with 35 Ncm (except for CreativeLine). CreativeLine should be fixed with 15 Ncm. The abutment screws may not be cemented or glued within the implant.
- Please observe the material specific minimum wall thickness for processing of the prosthetic abutments. The wall thickness may only be reduced partially – for titanium to 0.3 mm, for precious metal to 0.4 mm and for zirconium oxide to 0.5 mm.
- Color markings are provided for the abutments made of titanium (complete or partially) in the area of the FourByFour connection and the implant analogs according to the color guide system (green: Ø 3.40; yellow: Ø 3.80; red: Ø 4.10; blue: Ø 5.00).
- The connection surfaces of the abutments to the implant may not be blasted or treated. The hazard of an invasion of bacteria as well as weakening of the components must be avoided.

Model preparation: Impression copings / Implant analogs [Material Titanium grade 5]

Important: Impression Copings, Transfer Copings and Implant Analoges are for single use only. Multiple use is not allowed due to technical reasons (loss of precision of the thread parts resp. of Implant-Abutment-Connection). For preparation of a master model in the dental lab, the pertaining implant analogs are screw-fixed to the impression copings after impression taking. We distinguish between an open and a closed impression.

- **Open impression:** The impression copings are safely positioned and remained in the impression material of the (individually prepared) impression tray. Diameter-conform implant analogs are placed onto the FourByFour® Connection of the impression copings and screw-fixed with the long impression coping screw through the perforation of the individual tray.

- **Closed impression:** After curing of the impression material, the individually prepared impression tray is removed from the mouth of the patient. The impression copings are removed from the implants and screwed onto the corresponding implant analogs and then repositioned into the transfer copings which have remained in the impression. The transfer copings are for one-way use only. We recommend to prepare a removable gingival mask. The model should be made from high-quality class 4 plaster.

Temporary abutment “CreativeLine” [Material: Titanium grade 5]

Indications: Temporary crowns and bridges (time period is dependent on instructions of the acrylic material used but should not exceed 180 days)

Contraindications: Primary connections of abutments on diverging implants. General contraindications: see instructions for use of the implant system.

- The abutment is fixed with the abutment screw (FourByFour) safely on the implant analog.
- The knurled shaft of the abutment can be reduced in length according to the prosthetic situation.
- The circular groove marks the maximum length reduction.
- The modified abutment is then covered with opaque material.
- The slim funnel-shaped basic design favors a natural root profile in the surrounding gingival tissue by adding tooth-colored acrylic.
- The temporary abutment can be inserted immediately after exposure of the implant instead of a conventional non-rotation-safe titanium gingival former (healing abutment).
- A temporary crown can be inserted if the abutment is designed appropriately.

Note: If the CreativeLine abutment is inserted during the healing phase of the implant, any loading of the implant should be kept at a minimum. This is possible by connection of the temporary bridge with neighboring teeth, or by a corresponding reduction of the anatomy for avoiding load on the implants during occlusion and articulation.

Titanium abutment „VersaLine“ [Material: Titanium grade 5]

Indications: Preparation of primary parts in the telescope or conical crown technique, Cemented crown and bridge restorations

Contraindications: Primary connections of abutments, Prosthetic angulation to implant axis of more than 25°, Cast-on technique, Double crowns on implant diameter 3.40 mm, Single tooth restorations with free-end cantilever

- It is possible to prepare individual angulations between 0° and nearly 25° or to prepare circular shoulders which adjust to the natural contour of the gingiva.
- Furthermore, it is possible to veneer the abutment directly with an adequate titanium ceramics after special modification.
- We recommend to use titanium milling instruments for the preparation of this abutment.
- Different abutment variations are available:
 - In height of 7.5 mm with wide subgingival shape for low available gingiva
 - In height of 11.0 mm with a slim subgingival shape for high available gingiva

Important note: Do not use implants with a diameter of 3.40 mm for the restoration of telescope or conical crowns.

Titanium abutment “NaturalLine” [Material: Titanium grade 5]

Indications: Cemented crown and bridge restorations, Lateral screw-connected crowns and bridges (for individually prepared screw connections)

Contraindications: Primary connections of abutments, Prosthetic angulation to the implant axis of more than 25°, Cast-on technique, Double crowns on implant diameter 3.40 mm OT-F²/OT-F³, Single tooth restoration with free-end cantilever

- The distinguished feature is the subgingival design. Starting from the implant shoulder, the short version of the abutment (gingival height GH 1.2 mm) it extends circularly convex and ends in a surrounding gently curved chamfer.
- The higher version (gingiva height 3.5 mm) extends circularly concave from the implant and then forms a convex shape supporting the surrounding gingiva. The surrounding shoulder descends from the oral side to the esthetic zone.
- The massive contour allows a special reduction by using suitable instruments (titanium drill, polisher) and to design an optimal emergence profile.
- The circular shoulder enables an exact transition to the crown to be prepared. The shoulder can be adjusted to the gingival contour.
- For single crowns an anti-rotation cementation should be made (especially for 0° Natural-Line abutments). It may be necessary to prepare an additional rotation protection into the abutment at an early stage.
- Different abutment variations are available:
 - In angulations 0°, 15° and 25°
 - In socket heights of 1.20 mm and 3.5 mm

Goldbase Abutment “GoldLine” [Material: Base made of Gold/Platinum alloy Ceramcor, manufacturer C&M, Switzerland. Furnace: acrylic (POM)]

Indications: Individual abutments for fixed crown and bridge restorations, Individual primary parts in the telescope and conical crown technique

Contraindications: Single tooth restorations with free-end cantilever, Prosthetic angulation to the implant axis of more than 25°, Primary connections of abutments

- The base of the abutment consists of a highly precise manufactured cast-on Gold-platinum alloy and shows a perfect fit to the implant.
- The furnace fixed on the base is made of residue-free combustible acrylic.

Processing: The acrylic furnace of the abutment is shortened according to the prosthetic situation. The selected shape of the abutment is waxed onto the acrylic furnace. When preparing bridge or double crown constructions, a mutual insertion direction prior to modelation should be determined. Take care to prepare a circular functional surface of minimum 5 mm when designing the primary parts for the double crown technique. When investing the model, please be careful to keep a clear flow of the investing compound (containing phosphate) free of bubbles within the screw channel. Please note the instructions of the manufacturer of the investment compound. The final temperature for preheating should be kept for approx. 30 – 60 minutes according to the size of the muffle. A precious metal alloy is cast onto the existing precious metal basis using a regular process. The precious metal alloy should have a high level of consistency and observe the instructions of the manufacturer of the alloy. Leave the muffle to cool off at room temperature. In order to minimize the mechanical stress when removing the invested material, the precious metal base should be carefully polished with blast pearls at a pressure of approx. 1.5 bar. We recommend to use acid for removal of remaining investing material especially in the area of the implant-abutment interface (FourByFour), or to steam-clean or to remove in an ultrasonic bath. During further processing of the GoldLine abutment (drilling, blasting, polishing) the connection interface must be protected from damage, therefore the abutment should be screwed onto an implant analog. When preparing a ceramics crown for direct veneering with occlusal screw connection care must be taken that a minimum wall thickness of 0.5 mm is remaining. Due to the different thermal expansion coefficient (TEC) of the cast-on GoldLine base and the veneer ceramics, these should not have any contact. Please observe the instructions of the veneer ceramics manufacturer.

Important Note: Do not use implants of diameter 3.40 mm for telescopic or conical crown restorations !

Note for preparation:
Melting range: 1400 – 1490°C (2550-2710°F)
Cast-on temperature: up to 1350°C (2,462°F)
Thermal expansion coefficient (25-500°C): 11.9x10⁻⁶ K⁻¹,
Thermal expansion coefficient (25-600°C): 12.2x10⁻⁶ K⁻¹

Zirconium abutment “CeraLine” [Material: Base Titanium grade 5, Abutment Zirconium oxide]

Indications: Cemented crown and bridge restorations in full ceramics, Preparation of primary parts in the telescopic and conical crown technique

Contraindications: Primary connections of abutments, Prosthetic angulation to the implant axis of more than 25°, Cast-on technique, Double crowns on implant diameter 3.40 mm, Single tooth restoration with free-end cantilever

- The CeraLine is an abutment series of zirconium oxide abutment moulds which can be individualized.
- The connection to the implant is guaranteed by a highly precise titanium base. The abutment screw transmits the tear strength to the titanium base during final fixation, and not onto the zirconium part of the final abutment.
- The abutment is made of highly solid Yttrium stabilized zirconium oxide. These prefabricated moulds can be modified by trimming with a water cooled turbine.
- A veneering with suitable zirconium ceramics and color characterization is also possible.
- After these modifications have been made, the final zirconium part is fixed onto the titanium base by adhesive. The design of the base allows a precise positioning at adhesion so that transmission failures are avoided.

Processing: At first, the zirconium part is modified according to the prosthetic situation. Please use only water-cooled turbine handpieces and suitable diamond trimmers. When preparing a bridge or double crown construction, a mutual insertion direction should be determined first. A minimum wall thickness of 0.5 mm should remain at the final trimming of the zirconium part.

Adhesion: Basically the surfaces of both parts for adhesion should first be blasted carefully with 50µm aluminum oxide. Protect the FourByFour connection of the titanium base from damage by screwing the basis onto an implant analog prior to blasting and adhesion. Use only the laboratory screw and close the screw entrance of the titanium base with wax prior to adhesion. For adhesion please observe the instructions for use of the adhesive by the manufacturer.

- Suitable adhesives are among others:**
- Panavia F 2.0 OP (Opaque); Kuraray (permission of the manufacturer for using the Alloy Primer and the Ceramic Primer)
 - Relyx™ Unicem white Opaque; 3M Espe (permission of the manufacturer so far only for Lava™ Zirconium oxide)
 - Multilink Implant / Multilink Hybrid Abutment; Ivoclar Vivadent

Please make sure to perform an additional light curing.

Bar connection “ProfiLine” [Material: Adapter Titanium grade 5, Bar connector Titanium grade 5, Bar connector precious metal (Ceramcor®)]

Indications: Confectioned or individually milled bar constructions for anchoring of implant-supported complete prosthetic constructions in the maxilla and mandible.

Contraindications: Combination of periodontal (tooth) and implant-supported constructions, Combinations with other anchoring elements (for instance ball head attachment)

- The ProfiLine abutment is two-piece and serves to prepare confectioned as well as individually milled bar constructions.
- 1.) First step is to select the adapter corresponding to the implant diameter and the prevailing gingival height. The feature of this adapter is the internal contour being identical with the implant. This allows different procedures of preparation.
- This adapter can be screwed directly onto the implant in the mouth of the patient immediately by using a 4X4 insertion driver and fixed with a wrench at 35 Ncm. Since the internal contours are identical, the impression can be taken from adapter level, then provided with a healing abutment which remains within the patient’s mouth.
- The selection of the adapter is made by the dental technician and screwed onto the model implant analog in the laboratory.
- If the gingival height is small, the adapter is not used, and the bar connector can be placed directly onto the implant.
- 2.) Selecting the bar connector:
 - According to the type of the next following processing (cast-on, soldering, laser welding, adhesion) the bar connector either in titanium or precious metal should be used.
 - A compensation of divergency of the implants up to 26° is possible due to the conical connection between adapter and bar connector. If there are larger divergencies, we recommend to use a separating attachment.

- Note:**
- Please observe to keep an absolutely tension-free fit for all bar constructions (passive fit; Sheffield test) in the mouth of the patient.
 - Bilateral extensions of the bar beyond the posterior implants may be useful for static reasons. These should not, however, extend over a length of max. 5 mm (premolar width). Please consider the number of implants inserted for such extensions.
 - If only two implants are placed, the bar shape (round bar; drop-shaped) should permit a slight rotation of the anchored prosthetic construction parallel to the jaw articulation axis.
 - An off-central loading of the implants by an unfavorable bar connection is to be avoided.
 - Please observe basically to design a suitable dimension of the bar construction.
 - When using an adhesive, please follow the instructions of the manufacturer of the adhesive.

Bar Connector Titanium (grade 5)

- To be used for gluing or for laser welding with the individual casted titanium bar.
- Individually confectioned bars have to be separated from the bar connectors prior to investing and then casted separately. We recommend a welding on the master model or better an adhesion in the mouth of the patient. Please observe the instructions of the manufacturer of laser units. The bar connectors of titanium are not suitable for a direct cast-on.

Bar Connector Precious Metal (Ceramcor®)

- Confectioned bars can be directly soldered. Individually prepared bars can be casted on directly.
- We recommend to invest and to cast the bar modellation separately from the precious metal elements. This enables the user to solder, weld or glue the parts.
- The investment of the modellation is made with a phosphate-containing investment material. The instructions of the manufacturer of the investment material must be followed. The final temperature for preheating should be held for 30 – 60 minutes according to the size of the muffle.

Cast-on: In general processes the precious metal alloy can be directly casted onto the prefabricated precious metal bar connector. When selecting the precious metal alloy the highest

fixation abilities should be considered. Please follow the instructions of the alloy manufacturer. The muffle should slowly cool off to room temperature. In order to minimize mechanical load at removal from investment, the area of the precious metal bar connector should be carefully blasted with beads by using a pressure of 5.1.5 bar. We recommend to either use acid or steam or an ultrasonic bath for removing any remaining investment material especially in the area of the connection to the ProfiLine Adapter (cone). For the following preparation of the ProfiLine abutment (milling, blasting, polishing) the connection area should be protected from possible damage by placing it on a model analog.

Note for preparation:
Melting range: 1400 – 1490°C (2550-2710°F)
Cast-on temperature: up to 1350°C (2,462°F)

Adhesion:
We recommend to separate individually prepared bars from the bar connectors prior to investing and to cast them separately. For achieving a passive fit, we recommend to weld/solder on the master model or rather use an adhesive in the mouth of the patient.

Ball Head Abutment “TecLine” [Material: Ball abutment titanium grade 5, O-ring housing titanium grade 4, O-ring red: VMQ 70 (Methyl-vinyl-silicone) FKM 75, matrix Dalbo® Plus elliptic titanium grade 4, lamellae retention insert Eilitor®]

Indications: Anchorage of full prostheses in the edentulous maxilla and mandible in connection with two implants for achieving a transversal rotation axis, Anchorage of total prostheses in the edentulous maxilla and mandible in connection with four or also with six implants

Contraindications: An uneven number of implants per jaw, Implants inserted in non-symmetrical position, A combination with other retention elements (for instance telescopes or also natural teeth), Disparallel implants of more than 20° divergence

- This abutment is screwed onto the implant with the octagon driver.
- For anchoring, the user can choose between O-Ring attachment and Retention Anchor Dalbo® Plus elliptic*.

O-Ring Attachment

- At first a red O-Ring is inserted into a titanium metal housing. This ring remains in the housing during the laboratory processing and is replaced by the second red O-Ring for final insertion of the denture. As an option a black O-Ring is available for higher retention.

Retention Anchor Dalbo® Plus elliptic*

- This retention cap consists of two parts: a titanium housing with retention wings for fixation in the denture, and therein a Lamellae Retention Insert screw-in, made of precious metal (Eilitor®*), for which the forces can be individually adjusted by using an Activator Key (200-1,200 gram). Dalbo® Plus elliptic* can be used for up to 20° divergence of the implant. For detailed information, please refer to the separate instructions supplied with the product.

- Note:**
- **The implants should be positioned in the corresponding jaw area symmetrically, for static reasons and for achieving a parallel rotation axis to the jaw articulation axis. An uneven number of implants should be avoided.**
 - **Please select the smallest possible gingival height (GH) in order to avoid possible leverage load to the construction.**

* Manufacturer: Cendres & Métaux; Switzerland

Explanation of the symbols	
	Do not reuse
	Use by
	Keep dry
	Manufacturer
	Caution
	Do not use if the packaging is damaged
	Consult instruction for use
	None sterile